

astik's

Cryotube with screw cap Basic Line

Made of polypropylene with silicone-rubber lined cap. Conical bottom and skirt. With exterior screw. Can withstand temperatures from -70 °C up to +120 °C



reference	capacity	∅ (mm)	∅ stopper (mm)	h (mm)	type	pcs/pack	price/pack	average price/pc
CRTP-018-500	1,8 ml	11	13	46	graduated, blue cap	500		
CRTP-R18-500	1,8 ml	11	13	46	graduated, red cap	500		
CRTP-040-500	4 ml	12,5	15	53	without graduation	500		
"Saving package"								
CRTP-018-10K	1,8 ml	11	13	46	graduated, blue cap	20 x 500		
CRTP-R18-10K	1,8 ml	11	13	46	graduated, red cap	20 x 500		

astik's

Cryotube with screw cap, sterile, Standard Line

PP tube/HDPE Closure. Medical grade resin compliant with USP class VI. Leak proof with easy 3/4 turn and final 1/4 turn for perfect seal. Packed in resealable bags. Highly visible graduations and labeling area. Star foot design. Certified DNase, RNase and pyrogen-free. They withstand temperatures from -196 °C up to 40 °C. With exterior screw. Sterilized by gamma irradiation



reference	capacity	∅ (mm)	h (mm)	type	pcs/pack	price/pack	average price/pc
CRTT-010-025	1,0 ml	12	35	printed graduation, sterile	25		
CRTT-018-025	1,8 ml	12	48	printed graduation, sterile	25		
CRTT-045-010	4,5 ml	12	90	printed graduation, sterile	10		

astik's

Cryotube with screw cap, sterile, Premium Line

Both the lid and the tube are autoclavable; a small o-ring inside the lid ensures that the tube is tightly sealed (o-rings made from other materials may pose a risk of contamination in long-term storage). They withstand temperatures from -196 °C up to +120 °C. Certified DNase, RNase and endotoxin-free. With exterior screw. With graduation marks, a writing surface and an unique barcode printed directly on the tube. Their particular design allows them to be selfstanding. Sterilized by beta irradiation



reference	capacity	h (mm)	type	pcs/pack	price/pack	average price/pc
CRTD-018-050	1,8 ml	44,5	printed graduation, sterile	50		

astik's

Cryotubes with inner thread for 10x10 cryoboxes, Premium Line

Straightline design allows to be stored in standard 10x10 cryo boxes giving 23% additional storage capacity. Single-handed quick seal, closure allows tightening or removal with just half turn avoiding long-term injuries to the user. Star food stand makes it compatible with common working stations on the market. Specially designed cap with lip seal feature, 100% leak proof. Cap and tube are made of FDA-approved medical grade polypropylene with the same coefficient of expansion enhancing leak proof qualities within the temperature range (-196° to +120°C). Autoclavable. RNase and DNase free. Manufactured in a clean room environment, ISO-9001, ISO-14001 and ISO-13485. Sterilized according to ISO 11137-2:2013 "Sterility Assurance Level" (SAL) by standard products 10⁻⁶



reference	capacity	∅ (mm)	h (mm)	type	pcs/pack	price/pack	average price/pc
CRTH-N05-100	0,5 ml	12,08	28,6	non sterile, with writing area	100		
CRTH-N10-100	1 ml	12,08	36,1	non sterile, with writing area	100		
CRTH-N20-100	2 ml	12,08	47,6	non sterile, with writing area	100		
CRTH-S20-100	2 ml	12,08	47,6	sterile, with writing area	100		
CRTH-N40-100	4 ml	12,08	73,3	non sterile, with writing area	100		
CRTH-N50-100	5 ml	12,08	88,5	non sterile, with writing area	100		
"Saving package"							
CRTH-N20-1K0	2 ml	12,08	47,6	non sterile, with writing area	10 x100		
CRTH-S20-1K0	2 ml	12,08	47,6	sterile, with writing area	10 x100		

NEW

**Pack of coloured inserts for cryotubes with inner thread**

For Premium Line cryotubes, reference CRTH. 100 inserts of each colour

reference	tipo	pcs/pack	price/pack	average price/pc
CRTI-MCP-500	red, blue, green, yellow and purple	500		

Carton freezer box for cryotubes up to 2 ml

Delivered already assembled

reference	h (mm)	l (mm)	w (mm)	type	pcs/pack	price/pack
CBOX-081-001	47	134	134	up to 81 tubes ($\varnothing_{\text{max.}}$ 14 mm) white	1	
CBOX-R81-001	47	134	134	up to 81 tubes ($\varnothing_{\text{max.}}$ 14 mm) red	1	
CBOX-B81-001	47	134	134	up to 81 tubes ($\varnothing_{\text{max.}}$ 14 mm) blue	1	
CBOX-G81-001	47	134	134	up to 81 tubes ($\varnothing_{\text{max.}}$ 14 mm) green	1	
CBOX-Y81-001	47	134	134	up to 81 tubes ($\varnothing_{\text{max.}}$ 14 mm) yellow	1	
CBOX-N81-001	47	134	134	up to 81 tubes ($\varnothing_{\text{max.}}$ 14 mm) orange	1	
CBOX-100-001	47	134	134	up to 100 tubes ($\varnothing_{\text{max.}}$ 12 mm) white	1	
CBOX-10R-001	47	134	134	up to 100 tubes ($\varnothing_{\text{max.}}$ 12 mm) red	1	
CBOX-10B-001	47	134	134	up to 100 tubes ($\varnothing_{\text{max.}}$ 12 mm) blue	1	
CBOX-10G-001	47	134	134	up to 100 tubes ($\varnothing_{\text{max.}}$ 12 mm) green	1	
CBOX-10Y-001	47	134	134	up to 100 tubes ($\varnothing_{\text{max.}}$ 12 mm) yellow	1	
CBOX-10N-001	47	134	134	up to 100 tubes ($\varnothing_{\text{max.}}$ 12 mm) orange	1	
"Saving package"						
CBOX-081-096	47	134	134	up to 81 tubes ($\varnothing_{\text{max.}}$ 14 mm) white	96	
CBOX-100-096	47	134	134	up to 100 tubes ($\varnothing_{\text{max.}}$ 12 mm) white	96	

Carton freezer box for 4 ml cryotubes

Delivered already assembled

reference	h (mm)	l (mm)	w (mm)	type	pcs/pack	price/pack
CBOX-082-001	76	134	134	up to 81 tubes ($\varnothing_{\text{max.}}$ 14 mm) white	1	
CBOX-101-001	76	134	134	up to 100 tubes ($\varnothing_{\text{max.}}$ 12 mm) white	1	

Plastic freezer box for cryotubes

Made of polypropylene. Autoclavable

reference	h (mm)	l (mm)	w (mm)	type	pcs/pack	price/pack
FBXP-081-001	51	128	128	for 81 tubes up to 1.8 ml ($\varnothing_{\text{max.}}$ 12 mm)	1	
FBXP-100-001	53	140	140	for 100 tubes up to 1.8 ml ($\varnothing_{\text{max.}}$ 12 mm)	1	

**Plastic freezer box for cryotubes with numerical grid on lid**


Made of polypropylene with alphanumeric numeration printed on lid and bottom

reference	h (mm)	l (mm)	w (mm)	type	pcs/pack	price/pack
FBXP-B02-001	50	130	130	blue, for 81 tubes up to 2 ml ($\varnothing_{\text{max.}}$ 12 mm)	1	
FBXP-G02-001	50	130	130	green, for 81 tubes up to 2 ml ($\varnothing_{\text{max.}}$ 12 mm)	1	
FBXP-R02-001	50	130	130	red, for 81 tubes up to 2 ml ($\varnothing_{\text{max.}}$ 12 mm)	1	
FBXP-Y02-001	50	130	130	yellow, for 81 tubes up to 2 ml ($\varnothing_{\text{max.}}$ 12 mm)	1	
FBXP-B10-001	50	130	130	blue, for 100 tubes up to 2 ml ($\varnothing_{\text{max.}}$ 12 mm)	1	
FBXP-G10-001	50	130	130	green, for 100 tubes up to 2 ml ($\varnothing_{\text{max.}}$ 12 mm)	1	
FBXP-R10-001	50	130	130	red, for 100 tubes up to 2 ml ($\varnothing_{\text{max.}}$ 12 mm)	1	
FBXP-Y10-001	50	130	130	yellow, for 100 tubes up to 2 ml ($\varnothing_{\text{max.}}$ 12 mm)	1	


astik's




bb



bb




astik's



astik's



astik's




Plastic freezer box for cryotubes with snap-on lid

Made of polypropylene. Autoclavable and suitable for freezing (withstand temperatures from -130 °C up to +130 °C). Box with 81 alphanumeric printed positions and suitable for cryotubes up to 2 ml, microtubes 1.5 or 2 ml and chromatography vials 12 x 32 mm. Lid fixed with the box with a Snap-on device. Available in two versions: magenta colour box and natural lid; or black box and lid for photosensitive products

reference	h (mm)	l (mm)	w (mm)	type	pcs/pack	price/pack
FBXP-K01-001	62	140	140	clear/magenta, for 81 tubes ($\varnothing_{\text{max.}}$ 12 mm)	1	
FBXP-K02-001	62	140	140	black/black, for 81 tubes ($\varnothing_{\text{max.}}$ 12 mm)	1	

astik's



PC freezer storage boxes for cryotubes with numeric graduation on lid

Made of Polycarbonate (PC). Autoclavable. Orange color. Withstands temperatures from -196 °C up to 120 °C. Printed grid on the lid


reference	h (mm)	l (mm)	w (mm)	type	pcs/pack	price/pack
FBPC-N02-001	52	133	133	for 81 tubes up to 2 ml ($\varnothing_{\text{max.}}$ 12 mm)	1	
FBPC-N04-001	82	133	133	for 81 tubes up to 4 ml ($\varnothing_{\text{max.}}$ 12 mm)	1	
FBPC-N05-001	95	133	133	for 81 tubes up to 5 ml ($\varnothing_{\text{max.}}$ 12 mm)	1	

clink

Racks for upright freezers

Made of stainless steel. For boxes of 135 x 135 mm. Delivery time 2 weeks

reference	h (mm)	l (mm)	w (mm)	type	pcs/pack	price/pack
RACV-33A-001	165	420	140	for 3 x 3 cryoboxes of 50 mm tall	1	
RACV-43A-001	223	420	140	for 4 x 3 cryoboxes of 50 mm tall	1	
RACV-53A-001	277	420	140	for 5 x 3 cryoboxes of 50 mm tall	1	
RACV-44A-001	223	560	140	for 4 x 4 cryoboxes of 50 mm tall	1	
RACV-54A-001	277	560	140	for 5 x 4 cryoboxes of 50 mm tall	1	
RACV-33B-001	245	420	140	for 3 x 3 cryoboxes of 75 mm tall	1	
RACV-34B-001	245	560	140	for 3 x 4 cryoboxes of 75 mm tall	1	
RACV-23C-001	200	420	140	for 2 x 3 cryoboxes of 95 mm tall	1	
RACV-24C-001	200	560	140	for 2 x 4 cryoboxes of 95 mm tall	1	
RACV-23D-001	265	420	140	for 2 x 3 cryoboxes of 125 mm tall	1	
RACV-24D-001	265	560	140	for 2 x 4 cryoboxes of 125 mm tall	1	



clink

Racks for upright freezers with removable trays

Made of stainless steel. For boxes of 135 x 135 mm. The removable trays allow to use the carton boxes without having to remove the complete rack. Delivery time 2 weeks

reference	h (mm)	l (mm)	w (mm)	type	pcs/pack	price/pack
RACE-33A-001	170	440	143	for 3 x 3 cryoboxes of 50 mm tall	1	
RACE-43A-001	230	440	143	for 4 x 3 cryoboxes of 50 mm tall	1	
RACE-53A-001	285	440	143	for 5 x 3 cryoboxes of 50 mm tall	1	
RACE-44A-001	230	580	143	for 4 x 4 cryoboxes of 50 mm tall	1	
RACE-54A-001	285	580	143	for 5 x 4 cryoboxes of 50 mm tall	1	
RACE-33B-001	250	440	143	for 3 x 3 cryoboxes of 75 mm tall	1	
RACE-34B-001	250	580	143	for 3 x 4 cryoboxes of 75 mm tall	1	
RACE-23C-001	210	440	143	for 2 x 3 cryoboxes of 95 mm tall	1	
RACE-24C-001	210	580	143	for 2 x 4 cryoboxes of 95 mm tall	1	
RACE-23D-001	275	440	143	for 2 x 3 cryoboxes of 125 mm tall	1	
RACE-24D-001	275	580	143	for 2 x 4 cryoboxes of 125 mm tall	1	



Racks for chest freezers

Made of stainless steel. For boxes of 135 x 135 mm. Delivery time 2 weeks

reference	h (mm)	l (mm)	w (mm)	type	pcs/pack	price/pack
RACH-05A-001	280	140	140	for 5 cryoboxes of 50 mm tall	1	
RACH-07A-001	390	140	140	for 7 cryoboxes of 50 mm tall	1	
RACH-08A-001	445	140	140	for 8 cryoboxes of 50 mm tall	1	
RACH-09A-001	500	140	140	for 9 cryoboxes of 50 mm tall	1	
RACH-10A-001	555	140	140	for 10 cryoboxes of 50 mm tall	1	
RACH-11A-001	610	140	140	for 11 cryoboxes of 50 mm tall	1	
RACH-12A-001	665	140	140	for 12 cryoboxes of 50 mm tall	1	
RACH-13A-001	720	140	140	for 13 cryoboxes of 50 mm tall	1	
RACH-05B-001	405	140	140	for 5 cryoboxes of 75 mm tall	1	
RACH-06B-001	485	140	140	for 6 cryoboxes of 75 mm tall	1	
RACH-07B-001	565	140	140	for 7 cryoboxes of 75 mm tall	1	
RACH-08B-001	645	140	140	for 8 cryoboxes of 75 mm tall	1	
RACH-09B-001	725	140	140	for 9 cryoboxes of 75 mm tall	1	
RACH-04C-001	400	140	140	for 4 cryoboxes of 95 mm tall	1	
RACH-05C-001	500	140	140	for 5 cryoboxes of 95 mm tall	1	
RACH-06C-001	600	140	140	for 6 cryoboxes of 95 mm tall	1	
RACH-07C-001	700	140	140	for 7 cryoboxes of 95 mm tall	1	
RACH-04D-001	525	140	140	for 4 cryoboxes of 125 mm tall	1	
RACH-05D-001	655	140	140	for 5 cryoboxes of 125 mm tall	1	

clink



We can supply tailor-made sized racks for freezers. Kindly contact our Customer Service.

“Crushed Ice Maker

Compact crushed ice, perfectly suitable for the chemical and pharmaceutical sectors. Exterior in AISI 304 stainless steel and monoblock body in stainless steel with CFC's-free, high-pressure injected polyurethane insulation. Tank and door in ABS food-grade plastic. Production in optimum conditions: room temperature 21°C, water temperature 15°C. Water inlet Ø3/4 and outlet Ø24 mm. Cooling gas R404A. 220-240V/50Hz. 2 years warranty Delivery 4 weeks

reference	external dimensions (l x w x h)	production capacity	storage capacity	pcs/pack	price/pack
ICEM-060-001	450 x 620 x 680 mm	60 kg/24 h	10 kg	1	
ICEM-090-001	500 x 660 x 800 mm	90 kg/24 h	30 kg	1	

InfriCo
medicare



Ice bucket with lid, polyurethane

Blue color. Durable and rigid. Unbreakable in normal use. Convenient handle makes it easily portable. Excellent insulating properties suitable for use with liquid nitrogen, dry ice, water ice, salt solution and acetone

reference	capacity	dimensions (mm)	pcs/pack	price/pack
ICEP-B01-001	2,5 L	283 x 283 x 130	1	
ICEP-B02-001	4,5 L	280 x 280 x 185	1	

astik's





“Lab Care refrigerator +2 / + 8°C

Lab Care refrigerator available in solid and glass door type. Product temperature stability range of approximately ± 1 °C and temperature uniformity less than ± 2 °C. Great size display, digital controller with precision of 0.1 °C. 3 types of acoustic alarms (temperature rise or fall, probe failure and open door), remote alarm signal - potential free contact. They are equipped as standard with removable shelves in Epoxy coated galvanised steel and adjustable in height (3 shelves of 460 x 400 mm for the model of 150L and 5 shelves of 490 x 400 mm for the model of 360L), LED lighting and locks. Exterior / interior finishing: white plastic coated galvanised steel. They have multiple optional functions such as Bluetooth Datalogger, freezing protection controller, battery or RS485 communication port Modbus with PC or network, etc. 230V / 50Hz. 2 years warranty Delivery 4 weeks

reference	capacity	external dimensions (l x w x h)	shelves n°	door type	pcs/pack	price/pack
LCRE-15S-001	150 L	600 x 600 x 825 mm	3	solid	1	
LCRE-36S-001	360 L	620 x 600 x 1845 mm	5	solid	1	
LCRE-36G-001	360 L	620 x 600 x 1845 mm	5	glass	1	



“Lab Care Plus Freezer - 10 / -25°C

Lab Care freezer available in solid door type. Product temperature stability range of approximately ± 1 °C and temperature uniformity less than ± 2 °C. Great size display, digital controller with precision of 0.1 °C. 3 types of acoustic alarms (temperature rise or fall, probe failure and open door), remote alarm signal - potential free contact. They are equipped as standard with removable shelves in white plastic coated galvanised steel that allows you to adjust their position except for the model of 430 L that has 5 static chrome-plated wire shelf of 480 x 450 mm, LED lighting and locks. Exterior / interior finishing: white plastic coated galvanised steel. They have multiple optional functions such as Bluetooth Datalogger, battery or RS485 communication port Modbus with PC or network, etc. 230V / 50Hz. 2 years warranty Delivery 4 weeks

reference	capacity	external dimensions (l x w x h)	shelves n°	pcs/pack	price/pack
LCFR-15S-001	150 L	600 x 600 x 825 mm	3	1	
LCFR-40S-001	400 L	670 x 660 x 1955 mm	8	1	
LCFR-43S-001	430 L	687 x 670 x 2010 mm	5	1	



“Optional equipment for Lab Care refrigerators and Lab Care Plus Freezers

It is possible to customize the equipment depending on your needs and preferences

reference	description	pcs/pack	price/pack
LCAC-001-001	shelves for LCRE-15S, dim: 460 x 400 mm	1	
LCAC-002-001	shelves for LCRE-36S and -36G, dim: 490 x 400 mm	1	
LCAC-003-001	shelves for LCFR-15S, dim: 460 x 400 mm	1	
LCAC-004-001	shelves for LCFR-40S, dim: 455 x 425 mm	1	
LCAC-005-001	drawers for LCRE-15S, dim: 460 x 400 mm	1	
LCAC-006-001	drawers for LCRE-36S y -36G, dim: 490 x 400 mm	1	
LCAC-007-001	drawers for LCFR-15S, dim: 460 x 400 mm	1	
LCAC-008-001	drawers for LCFR-40S, dim: 455 x 425 mm	1	
LCAC-009-001	freezing protection controller for LCRE-15S, -36S y -36G	1	
LCAC-010-001	access port for all models	1	
LCAC-011-004	casters for all models	4	
LCAC-012-001	product simulation temperature probe for all models	1	
LCAC-013-001	battery for all models	1	
LCAC-014-001	RS485 communication port Modbus for all models	1	
LCAC-015-001	bluetooth datalogger for all models	1	

