

# **MICROORGANISMS**

## **LYOPHILIZED REFERENCE ORGANISMS**

**MICROSWAB**

**MICROSWABPLUS**

**MICROPELLETTRT**

**MICROPELLETTRTPLUS**

**MICROPELLETTRT10**

**PRODUCT LIST  
2019**



## **MECCONTI offers five different product lines of lyophilized micro-organisms.**

### **MICROSWABS**

A *MicroSwab* consists of a lyophilized pellet of a single micro-organism strain inside a pen-like device. This pen-like device contains a sterile swab and rehydration fluid to facilitate the transfer of the micro-organisms directly to culture media. The pellet has no defined CFU assay but a minimum of 1000 CFU/pellet to ensure proper growth. Use *MicroSwabs* whenever a defined assay value is not needed. The *MicroSwabs* are up to **4 passages removed from reference strain**.

### **MICROSWABPLUS**

The *MicroSwabPlus* also consists of a lyophilized pellet of a single micro-organism strain inside a pen-like device. This pen-like device contains a sterile swab and rehydration fluid to facilitate the transfer of the micro-organisms directly to culture media. The pellet has no defined CFU assay but a minimum of 1000 CFU/pellet to ensure proper growth. Whereas *MicroSwabs* are up to 4 passages removed from reference strain the *MicroSwabPlus* are only **2 passages removed from reference strain**.

### **MICROPELLETRT**

A *MicroPelletRT* consists of a lyophilized pellet of a single micro-organism strain with a defined mean assay value inside a plastic vial and one separate vial with 1 ml rehydration fluid. Adding the reconstitution fluid to the pellet vial it forms a suspension with a defined assay value of 10-100 CFU/0.1 ml as mandated by USP and Ph.Eur. for growth promotion tests.

### **MICROPELLETRT10**

*MicroPelletRT10* is a bigger format of *MicroPelletRT* and is intended for customers who perform a lot of testing.

A *MicroPelletRT10* consists of a lyophilized pellet of a single micro-organism strain with a defined mean assay value inside a plastic vial and one separate vial with 10 ml rehydration fluid. Adding pellet to the reconstitution fluid vial it forms a suspension with a defined assay value of 10-100 CFU/0.1 ml as mandated by USP and Ph. Eur. for growth promotion tests.

### **MICROPELETSRTPLUS**

*MicroPelletsRTPlus* microorganisms can be used to perform various essential quality control functions including microbial detection and enumeration, equipment calibration, method validation, bioburden determination, antibacterial effectiveness and lethality testing. *MicroPelletRTPlus* is available in concentrations ranging from  $10^2$  to  $10^8$  CFU per pellet. Rehydration fluid for this product line is available in sets of 2 vials of 1 ml each as a separate product.

**All products included in the list of microorganisms products ( pages 4-26)  
are available to order in two packaging configurations:**

MicroSwabs & MicroSwabsPlus - as set of 2 or 10 swabs per pack

MicroPelletsRT, MicroPelletsRTPlus & MicroPelletsRT10 - as set of 2 or 10 pellets per pack

# MICROSWAB

code	product name
<b>MSA001</b>	MicroSwab Acinetobacter baumannii ATCC® 19606™
<b>MSA002</b>	MicroSwab Acinetobacter Iwoffii ATCC® 17925™
<b>MSA003</b>	MicroSwab Acinetobacter sp. AmMS 203 ATCC® 49137™
<b>MSA004</b>	MicroSwab Acinetobacter sp. AmMs 202 ATCC® 49139™
<b>MSA005</b>	MicroSwab Acinetobacter sp. AmMS 243 ATCC® 49466™
<b>MSA006</b>	MicroSwab Actinomyces odontolyticus ATCC® 17929™
<b>MSA007</b>	MicroSwab Actinomyces pyogenes ATCC® 19411™
<b>MSA008</b>	MicroSwab Actinomyces pyogenes ATCC® 49698™
<b>MSA009</b>	MicroSwab Aerococcus viridans ATCC® 700406™
<b>MSA010</b>	MicroSwab Aeromonas hydrophila ATCC® 35654™
<b>MSA011</b>	MicroSwab Aeromonas hydrophila AmMS 199 ATCC® 49140™
<b>MSA012</b>	MicroSwab Aeromonas hydrophila WDCM 00063-ATCC® 7966™
<b>MSA013</b>	MicroSwab Aeromonas salmonicida ATCC® 33658™
<b>MSA014</b>	MicroSwab Alcaligenes faecalis ATCC® 35655™
<b>MSA015</b>	MicroSwab Alcaligenes faecalis ATCC® 8750™
<b>MSA016</b>	MicroSwab Aspergillus fumigatus ATCC® 204305™
<b>MSA017</b>	MicroSwab Aspergillus brasiliensis WDCM 00053-ATCC® 16404™
<b>MSA018</b>	MicroSwab Abiotrophia defectiva ATCC® 49176™
<b>MSA019</b>	MicroSwab Achromobacter xylosoxidans ATCC® 27061™
<b>MSA020</b>	MicroSwab Acinetobacter baumannii ATCC® BAA-747™
<b>MSA021</b>	MicroSwab Actinobacillus pleuropneumoniae ATCC® 27090™
<b>MSA022</b>	MicroSwab Actinomyces viscosus ATCC® 15987™
<b>MSB001</b>	MicroSwab Bacillus cereus ATCC® 10876™
<b>MSB002</b>	MicroSwab Bacillus cereus WDCM 00001-ATCC® 11778™
<b>MSB003</b>	MicroSwab Bacillus cereus ATCC® 14579™
<b>MSB004</b>	MicroSwab Paenibacillus polymyxa ATCC® 842™
<b>MSB005</b>	MicroSwab Lysinibacillus sphaericus ATCC® 4525™
<b>MSB006</b>	MicroSwab Geobacillus stearothermophilus ATCC® 7953™
<b>MSB007</b>	MicroSwab Bacillus subtilis subsp. spizizenii WDCM 00003-ATCC® 6633™
<b>MSB008</b>	MicroSwab Bacteroides distasonis ATCC® 8503™
<b>MSB009</b>	MicroSwab Bacteroides fragilis ATCC® 23745™
<b>MSB010</b>	MicroSwab Bacteroides fragilis ATCC® 25285™
<b>MSB011</b>	MicroSwab Bacteroides ovatus ATCC® 8483™
<b>MSB012</b>	MicroSwab Bacteroides thetaiomicron ATCC® 29741™



© The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalogue marks are trademarks of ATCC. Meccanti S.p.A.R.L. Sp. z o. o. is licensed to use these trademarks and sell products derived from ATCC cultures. Look for the ATCC Licensed Derivative Emblem for products derived from ATCC cultures.

code	product name
<b>MSB013</b>	MicroSwab Bacteroides uniformis ATCC® 8492™
<b>MSB014</b>	MicroSwab Bacteroides ureolyticus ATCC® 33387™
<b>MSB015</b>	MicroSwab Bacteroides vulgatus ATCC® 29327™
<b>MSB016</b>	MicroSwab Bacteroides vulgatus ATCC® 8482™
<b>MSB017</b>	MicroSwab Geotrichum capitatum ATCC® 10663™
<b>MSB018</b>	MicroSwab Bordetella bronchiseptica ATCC® 10580™
<b>MSB019</b>	MicroSwab Bordetella bronchiseptica ATCC® 4617™
<b>MSB020</b>	MicroSwab Bordetella pertussis ATCC® 8467™
<b>MSB021</b>	MicroSwab Moraxella catarrhalis ATCC® 23246™
<b>MSB022</b>	MicroSwab Moraxella catarrhalis ATCC® 25238™
<b>MSB023</b>	MicroSwab Moraxella catarrhalis ATCC® 25240™
<b>MSB024</b>	MicroSwab Moraxella catarrhalis ATCC® 49143™
<b>MSB025</b>	MicroSwab Moraxella catarrhalis ATCC® 8176™
<b>MSB026</b>	MicroSwab Brevundimonas diminuta ATCC® 19146™
<b>MSB027</b>	MicroSwab Burkholderia cepacia ATCC® 25416™
<b>MSB028</b>	MicroSwab Burkholderia cepacia ATCC® 25608™
<b>MSB029</b>	MicroSwab Ralstonia insidiosa ATCC® 49129™
<b>MSB030</b>	MicroSwab Bordetella parapertussis ATCC® 15311™
<b>MSB031</b>	MicroSwab Bacillus pumilus ATCC® 14884™
<b>MSB032</b>	MicroSwab Brochothrix thermosphacta WDCM 00071-ATCC® 11509™
<b>MSB034</b>	MicroSwab Bacteroides ovatus ATCC® BAA 1296™
<b>MSB035</b>	MicroSwab Bacillus circulans ATCC® 4513™
<b>MSB036</b>	MicroSwab Bacillus licheniformis WDCM 00068-ATCC® 14580™
<b>MSB037</b>	MicroSwab Bacillus megaterium ATCC® 6458™
<b>MSB038</b>	MicroSwab Bacillus pumilus ATCC® 21143™
<b>MSB039</b>	MicroSwab Bacillus amyloliquefaciens ATCC® 23842™
<b>MSB040</b>	MicroSwab Bacteroides ovatus ATCC® BAA1304™
<b>MSC001</b>	MicroSwab Campylobacter jejuni ATCC® 33292™
<b>MSC002</b>	MicroSwab Campylobacter jejuni WDCM 00156-ATCC® 29428™
<b>MSC003</b>	MicroSwab Campylobacter jejuni WDCM 00005-ATCC® 33291™
<b>MSC004</b>	MicroSwab Candida albicans ATCC® 36232™
<b>MSC005</b>	MicroSwab Candida albicans WDCM 00054-ATCC® 10231™
<b>MSC006</b>	MicroSwab Candida albicans ATCC® 14053™



The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalogue marks are trademarks of ATCC. Meccanti S.p.A. is licensed to use these trademarks and sell products derived from ATCC cultures. Look for the ATCC Licensed Derivative Emblem for products derived from ATCC cultures.

# MICROSWAB

code	product name
<b>MSC007</b>	MicroSwab Candida albicans ATCC® 60193™
<b>MSC008</b>	MicroSwab Candida albicans ATCC® 66027™
<b>MSC009</b>	MicroSwab Candida glabrata ATCC® 15126™
<b>MSC010</b>	MicroSwab Candida glabrata ATCC® 2001™
<b>MSC011</b>	MicroSwab Candida glabrata ATCC® 66032™
<b>MSC012</b>	MicroSwab Candida guilliermondii ATCC® 6260™
<b>MSC013</b>	MicroSwab Candida kefyr ATCC® 2512™
<b>MSC014</b>	MicroSwab Candida kefyr ATCC® 66028™
<b>MSC015</b>	MicroSwab Candida krusei ATCC® 14243™
<b>MSC016</b>	MicroSwab Issatchenka orientalis ATCC® 6258™
<b>MSC017</b>	MicroSwab Candida parapsilosis ATCC® 22019™
<b>MSC018</b>	MicroSwab Candida parapsilosis ATCC® 34136™
<b>MSC019</b>	MicroSwab Candida tropicalis ATCC® 66029™
<b>MSC020</b>	MicroSwab Candida tropicalis ATCC® 750™
<b>MSC021</b>	MicroSwab Elizabethkingia meningoseptica ATCC® 13253™
<b>MSC022</b>	MicroSwab Citrobacter freundii ATCC® 8090™
<b>MSC023</b>	MicroSwab Clostridium difficile ATCC® 9689™
<b>MSC024</b>	MicroSwab Clostridium histolyticum ATCC® 19401™
<b>MSC025</b>	MicroSwab Clostridium novyi ATCC® 7659™
<b>MSC026</b>	MicroSwab Clostridium perfringens WDCM 00007-ATCC® 13124™
<b>MSC028</b>	MicroSwab Clostridium sordellii ATCC® 9714™
<b>MSC029</b>	MicroSwab Clostridium sporogenes ATCC® 11437™
<b>MSC030</b>	MicroSwab Clostridium sporogenes WDCM 00008-ATCC® 19404™
<b>MSC031</b>	MicroSwab Clostridium sporogenes ATCC® 3584™
<b>MSC032</b>	MicroSwab Corynebacterium pseudodiphtheriticum ATCC® 10700™
<b>MSC033</b>	MicroSwab Corynebacterium pseudodiphtheriticum ATCC® 10701™
<b>MSC034</b>	MicroSwab Exiguobacterium aurantiacum ATCC® 49676™
<b>MSC035</b>	MicroSwab Corynebacterium xerosis ATCC® 373™
<b>MSC036</b>	MicroSwab Cryptococcus albidus ATCC® 10666™
<b>MSC037</b>	MicroSwab Cryptococcus albidus ATCC® 34140™
<b>MSC038</b>	MicroSwab Cryptococcus albidus ATCC® 66030™
<b>MSC039</b>	MicroSwab Cryptococcus laurentii ATCC® 66036™
<b>MSC040</b>	MicroSwab Cryptococcus laurentii ATCC® 76483™
<b>MSC041</b>	MicroSwab Cryptococcus luteolus ATCC® 10671™
<b>MSC042</b>	MicroSwab Cryptococcus neoformans ATCC® 14116™



The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalogue marks are trademarks of ATCC. Mecconti S.A.R.L. Sp. z o. o. is licensed to use these trademarks and sell products derived from ATCC cultures. Look for the ATCC Licensed Derivative Emblem for products derived from ATCC cultures.

# MICROSWAB

code	product name
<b>MSC043</b>	MicroSwab Cryptococcus neoformans ATCC® 32045™
<b>MSC044</b>	MicroSwab Cryptococcus neoformans ATCC® 34877™
<b>MSC045</b>	MicroSwab Cryptococcus neoformans ATCC® 66031™
<b>MSC046</b>	MicroSwab Cryptococcus neoformans ATCC® 76484™
<b>MSC047</b>	MicroSwab Cryptococcus uniguttulatus ATCC® 66033™
<b>MSC048</b>	MicroSwab Cryptococcus humicolus ATCC® 9949™
<b>MSC049</b>	MicroSwab Candida albicans ATCC® 90028™
<b>MSC050</b>	MicroSwab Campylobacter coli WDCM 00004-ATCC® 43478™
<b>MSC051</b>	MicroSwab Citrobacter freundii WDCM 00006-ATCC® 43864™
<b>MSC052</b>	MicroSwab Candida glabrata ATCC® 90030™
<b>MSC053</b>	MicroSwab Clostridium septicum ATCC® 12464™
<b>MSC054</b>	MicroSwab Candida albicans ATCC® 64548™
<b>MSC055</b>	MicroSwab Chryseobacterium indologenes ATCC® 49471™
<b>MSC056</b>	MicroSwab Citrobacter freundii ATCC® 8454™
<b>MSC057</b>	MicroSwab Clostridium ramosum ATCC® 13937™
<b>MSC058</b>	MicroSwab Campylobacter jejuni ATCC® 49943™
<b>MSC059</b>	MicroSwab Corynebacterium striatum ATCC® BAA1293™
<b>MSC060</b>	MicroSwab Cryptococcus laurentii ATCC® 18803™
<b>MSC061</b>	MicroSwab Cardiobacterium hominis ATCC® 14900™
<b>MSD001</b>	MicroSwab Delftia acidovorans ATCC® 15668™
<b>MSE001</b>	MicroSwab Enterobacter aerogenes WDCM 00175-ATCC® 13048™
<b>MSE002</b>	MicroSwab Enterobacter aerogenes ATCC® 35028™
<b>MSE003</b>	MicroSwab Enterobacter cloacae WDCM 00083-ATCC® 13047™
<b>MSE004</b>	MicroSwab Enterobacter cloacae WDCM 00082-ATCC® 23355™
<b>MSE005</b>	MicroSwab Enterobacter cloacae ATCC® 49141™
<b>MSE006</b>	MicroSwab Enterobacter gergoviae ATCC® 33028™
<b>MSE007</b>	MicroSwab Enterococcus avium ATCC® 14025™
<b>MSE008</b>	MicroSwab Enterococcus avium ATCC® 49464™
<b>MSE009</b>	MicroSwab Enterococcus durans ATCC® 11576™
<b>MSE010</b>	MicroSwab Enterococcus durans ATCC® 19432™
<b>MSE011</b>	MicroSwab Enterococcus durans ATCC® 49135™
<b>MSE012</b>	MicroSwab Enterococcus durans ATCC® 49479™
<b>MSE013</b>	MicroSwab Enterococcus durans ATCC® 6056™



The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalogue marks are trademarks of ATCC. Meccanti S.p.A. is licensed to use these trademarks and sell products derived from ATCC cultures. Look for the ATCC Licensed Derivative Emblem for products derived from ATCC cultures.

# MICROSWAB

code	product name
<b>MSE014</b>	MicroSwab Enterococcus faecalis WDCM 00009-ATCC® 19433™
<b>MSE015</b>	MicroSwab Enterococcus faecalis WDCM 00087-ATCC® 29212™
<b>MSE016</b>	MicroSwab Enterococcus faecalis ATCC® 49149™
<b>MSE017</b>	MicroSwab Enterococcus faecalis WDCM 00085,WDCM 00152-ATCC® 51299™
<b>MSE018</b>	MicroSwab Enterococcus faecium ATCC® 35667™
<b>MSE019</b>	MicroSwab Enterococcus gallinarum ATCC® 700425™
<b>MSE020</b>	MicroSwab Epidermophyton floccosum ATCC® 52066™
<b>MSE021</b>	MicroSwab Erysipelothrix rhusiopathiae ATCC® 19414™
<b>MSE022</b>	MicroSwab Escherichia coli ATCC® 11303™
<b>MSE023</b>	MicroSwab Escherichia coli WDCM 00090-ATCC® 11775™
<b>MSE024</b>	MicroSwab Escherichia coli ATCC® 14948™
<b>MSE025</b>	MicroSwab Escherichia coli WDCM 00013-ATCC® 25922™
<b>MSE027</b>	MicroSwab Escherichia coli ATCC® 35218™
<b>MSE028</b>	MicroSwab Escherichia coli ATCC® 35421™
<b>MSE030</b>	MicroSwab Escherichia coli ATCC® 51446™
<b>MSE031</b>	MicroSwab Escherichia coli ATCC® 51755™
<b>MSE032</b>	MicroSwab Escherichia coli WDCM 00012,WDCM 00196-ATCC® 8739™
<b>MSE033</b>	MicroSwab Eggerthella lenta ATCC® 43055™
<b>MSE034</b>	MicroSwab Cronobacter muytjensii WDCM 00213- ATCC® 51329™
<b>MSE035</b>	MicroSwab Enterococcus faecalis WDCM 00210-ATCC® 33186™
<b>MSE036</b>	MicroSwab Enterococcus hirae WDCM 00011-ATCC® 10541™
<b>MSE037</b>	MicroSwab Escherichia coli ATCC® 10536™
<b>MSE038</b>	MicroSwab Escherichia coli ATCC® 9637™
<b>MSE039</b>	MicroSwab Enterococcus casseliflavus ATCC® 700327™
<b>MSE040</b>	MicroSwab Enterococcus faecium WDCM 00177-ATCC® 6057™
<b>MSE046</b>	MicroSwab Eikenella corrodens ATCC® BAA1152™
<b>MSE047</b>	MicroSwab Escherichia coli ATCC® 11229™
<b>MSE048</b>	MicroSwab Enterobacter hormaechei ATCC® 700323™
<b>MSF001</b>	MicroSwab Myroides odoratus ATCC® 4651™
<b>MSF002</b>	MicroSwab Fusobacterium nucleatum ATCC® 25586™
<b>MSF003</b>	MicroSwab Fluoribacter bozemanae ATCC® 33217™
<b>MSG001</b>	MicroSwab Gardnerella vaginalis ATCC® 14018™



© The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalogue marks are trademarks of ATCC. Meconti S.p.A. is licensed to use these trademarks and sell products derived from ATCC cultures. Look for the ATCC Licensed Derivative Emblem for products derived from ATCC cultures.

code	product name
<b>MSG002</b>	MicroSwab Gardnerella vaginalis ATCC® 49145™
<b>MSG003</b>	MicroSwab Geobacillus stearothermophilus ATCC® 12976™
<b>MSG004</b>	MicroSwab Geobacillus stearothermophilus ATCC® 12980™
<b>MSH001</b>	MicroSwab Haemophilus aphrophilus ATCC® 19415™
<b>MSH002</b>	MicroSwab Haemophilus haemolyticus ATCC® 33390™
<b>MSH003</b>	MicroSwab Haemophilus influenzae ATCC® 10211™
<b>MSH004</b>	MicroSwab Haemophilus influenzae ATCC® 19418™
<b>MSH005</b>	MicroSwab Haemophilus influenzae ATCC® 33533™
<b>MSH006</b>	MicroSwab Haemophilus influenzae ATCC® 33930™
<b>MSH007</b>	MicroSwab Haemophilus influenzae ATCC® 35056™
<b>MSH008</b>	MicroSwab Haemophilus influenzae ATCC® 49144™
<b>MSH009</b>	MicroSwab Haemophilus influenzae ATCC® 49247™
<b>MSH010</b>	MicroSwab Haemophilus influenzae ATCC® 49766™
<b>MSH011</b>	MicroSwab Haemophilus influenzae ATCC® 9006™
<b>MSH013</b>	MicroSwab Haemophilus parahaemolyticus ATCC® 10014™
<b>MSH014</b>	MicroSwab Haemophilus parainfluenzae ATCC® 51505™
<b>MSH015</b>	MicroSwab Aggregatibacter aphrophilus ATCC® 7901™
<b>MSH016</b>	MicroSwab Haemophilus paraphrophilus ATCC® 49146™
<b>MSH017</b>	MicroSwab Haemophilus influenzae ATCC® 9833™
<b>MSH018</b>	MicroSwab Haemophilus influenzae ATCC® CL 511™
<b>MSH019</b>	MicroSwab Haemophilus paraphrophilus ATCC® 49917™
<b>MSK001</b>	MicroSwab Klebsiella oxytoca ATCC® 43086™
<b>MSK002</b>	MicroSwab Klebsiella oxytoca ATCC® 49131™
<b>MSK003</b>	MicroSwab Klebsiella pneumoniae ATCC® 10031™
<b>MSK004</b>	MicroSwab Klebsiella pneumoniae WDCM 00097-ATCC® 13883™
<b>MSK005</b>	MicroSwab Klebsiella pneumoniae ATCC® 23357™
<b>MSK006</b>	MicroSwab Klebsiella pneumoniae ATCC® 27736™
<b>MSK007</b>	MicroSwab Klebsiella pneumoniae ATCC® 33495™
<b>MSK008</b>	MicroSwab Klebsiella pneumoniae ATCC® 35657™
<b>MSK009</b>	MicroSwab Klebsiella pneumoniae ATCC® 700603™
<b>MSK010</b>	MicroSwab Klebsiella pneumoniae ATCC® 9997™
<b>MSK011</b>	MicroSwab Kocuria rosea ATCC® 186™
<b>MSK012</b>	MicroSwab Klebsiella oxytoca ATCC® 700324™



© The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalogue marks are trademarks of ATCC. Meccanti S.p.A. is licensed to use these trademarks and sell products derived from ATCC cultures. Look for the ATCC Licensed Derivative Emblem for products derived from ATCC cultures.

# MICROSWAB

code	product name
<b>MSK013</b>	MicroSwab Klebsiella pneumoniae ATCC® BAA1705™
<b>MSK014</b>	MicroSwab Klebsiella pneumoniae ATCC® BAA1706™
<b>MSK015</b>	MicroSwab Kingella kingae ATCC® 23330™
<b>MSK016</b>	MicroSwab Klebsiella variicola WDCM 00206-ATCC® 31488™
<b>MSK017</b>	MicroSwab Kocuria kristinae ATCC® BAA-752™
<b>MSK018</b>	MicroSwab Kurthia gibsonii ATCC® 43195™
<b>MSL001</b>	MicroSwab Lactobacillus acidophilus WDCM 00098-ATCC® 4356™
<b>MSL002</b>	MicroSwab Lactobacillus casei subsp. casei WDCM 00100-ATCC® 393™
<b>MSL003</b>	MicroSwab Lactobacillus rhamnosus WDCM 00101-ATCC® 7469™
<b>MSL004</b>	MicroSwab Lactococcus lactis subsp. lactis ATCC® 11454™
<b>MSL005</b>	MicroSwab Leclercia adecarboxylata ATCC® 23216™
<b>MSL006</b>	MicroSwab Legionella pneumophila WDCM 00107-ATCC® 33152™
<b>MSL007</b>	MicroSwab Listeria grayi ATCC® 25401™
<b>MSL008</b>	MicroSwab Listeria innocua WDCM 00017-ATCC® 33090™
<b>MSL010</b>	MicroSwab Listeria ivanovii WDCM 00018-ATCC® 19119™
<b>MSL011</b>	MicroSwab Listeria monocytogenes WDCM 00020-ATCC® 19111™
<b>MSL012</b>	MicroSwab Listeria monocytogenes ATCC® 7644™
<b>MSL013</b>	MicroSwab Listeria monocytogenes ATCC® 19115™
<b>MSL014</b>	MicroSwab Lactobacillus plantarum ATCC® 8014™
<b>MSL015</b>	MicroSwab Listeria monocytogenes ATCC® 15313™
<b>MSL023</b>	MicroSwab Listeria monocytogenes WDCM 00021-ATCC® 13932™
<b>MSL024</b>	MicroSwab Listeria monocytogenes ATCC® BAA751™
<b>MSL025</b>	MicroSwab Listeria monocytogenes ATCC® 19117™
<b>MSL026</b>	MicroSwab Listeria seeligeri ATCC® 35967™
<b>MSL027</b>	MicroSwab Listeria welshimeri ATCC® 35897™
<b>MSM001</b>	MicroSwab Micrococcus luteus WDCM 00111-ATCC® 4698™
<b>MSM002</b>	MicroSwab Micrococcus luteus ATCC® 49732™
<b>MSM003</b>	MicroSwab Kocuria rhizophila ATCC® 533™
<b>MSM004</b>	MicroSwab Kocuria rhizophila ATCC® 9341™
<b>MSM005</b>	MicroSwab Micrococcus sp. ATCC® 700405™
<b>MSM006</b>	MicroSwab Microsporum canis ATCC® 11621™
<b>MSM007</b>	MicroSwab Morganella morganii WDCM 00112-ATCC® 25830™
<b>MSM008</b>	MicroSwab Moraxella osloensis ATCC® 10973™
<b>MSM009</b>	MicroSwab Micrococcus luteus ATCC® 10240™



The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalogue marks are trademarks of ATCC. Mecconti S.p.A.R.L. Sp. z o. o. is licensed to use these trademarks and sell products derived from ATCC cultures. Look for the ATCC Licensed Derivative Emblem for products derived from ATCC cultures.

# MICROSWAB

code	product name
<b>MSM012</b>	MicroSwab Malassezia furfur ATCC® 14521™
<b>MSN001</b>	MicroSwab Neisseria gonorrhoeae ATCC® 19424™
<b>MSN002</b>	MicroSwab Neisseria gonorrhoeae ATCC® 31426™
<b>MSN003</b>	MicroSwab Neisseria gonorrhoeae ATCC® 43069™
<b>MSN004</b>	MicroSwab Neisseria gonorrhoeae ATCC® 43070™
<b>MSN005</b>	MicroSwab Neisseria gonorrhoeae ATCC® 49226™
<b>MSN006</b>	MicroSwab Neisseria gonorrhoeae ATCC® 49981™
<b>MSN007</b>	MicroSwab Neisseria lactamica ATCC® 23970™
<b>MSN008</b>	MicroSwab Neisseria lactamica ATCC® 23971™
<b>MSN009</b>	MicroSwab Neisseria lactamica ATCC® 49142™
<b>MSN010</b>	MicroSwab Neisseria meningitidis Group A ATCC® 13077™
<b>MSN011</b>	MicroSwab Neisseria meningitidis Group B ATCC® 13090™
<b>MSN012</b>	MicroSwab Neisseria meningitidis Group C ATCC® 13102™
<b>MSN013</b>	MicroSwab Neisseria mucosa ATCC® 19695™
<b>MSN014</b>	MicroSwab Neisseria sicca ATCC® 9913™
<b>MSN015</b>	MicroSwab Neisseria perflava ATCC® 14799™
<b>MSN016</b>	MicroSwab Nocardia asteroides ATCC® 19247™
<b>MSN017</b>	MicroSwab Nocardia brasiliensis ATCC® 19296™
<b>MSN018</b>	MicroSwab Nocardia farcinica ATCC® 3308™
<b>MSO001</b>	MicroSwab Oligella ureolytica ATCC® 43534™
<b>MSO002</b>	MicroSwab Oligella urethralis ATCC® 17960™
<b>MSO003</b>	MicroSwab Ochrobactrum anthropi ATCC® 49687™
<b>MSP001</b>	MicroSwab Pasteurella multocida ATCC® 43137™
<b>MSP002</b>	MicroSwab Penicillium chrysogenum ATCC® 10106™
<b>MSP003</b>	MicroSwab Peptostreptococcus anaerobius ATCC® 27337™
<b>MSP004</b>	MicroSwab Peptoniphilus asaccharolyticus ATCC® 29743™
<b>MSP005</b>	MicroSwab Finegoldia magna ATCC® 29328™
<b>MSP006</b>	MicroSwab Parvimonas micra ATCC® 33270™
<b>MSP007</b>	MicroSwab Plesiomonas shigelloides ATCC® 14029™
<b>MSP008</b>	MicroSwab Plesiomonas shigelloides ATCC® 51903™
<b>MSP009</b>	MicroSwab Porphyromonas gingivalis ATCC® 33277™
<b>MSP010</b>	MicroSwab Porphyromonas levii ATCC® 29147™
<b>MSP011</b>	MicroSwab Propionibacterium acnes ATCC® 11827™
<b>MSP012</b>	MicroSwab Proteus mirabilis ATCC® 12453™



The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalogue marks are trademarks of ATCC. Meccanti S.p.A. is licensed to use these trademarks and sell products derived from ATCC cultures. Look for the ATCC Licensed Derivative Emblem for products derived from ATCC cultures.

# MICROSWAB

code	product name
<b>MSP013</b>	MicroSwab <i>Proteus mirabilis</i> ATCC® 29245™
<b>MSP014</b>	MicroSwab <i>Proteus mirabilis</i> WDCM 00023-ATCC® 29906™
<b>MSP015</b>	MicroSwab <i>Proteus mirabilis</i> ATCC® 35659™
<b>MSP016</b>	MicroSwab <i>Proteus mirabilis</i> ATCC® 43071™
<b>MSP017</b>	MicroSwab <i>Proteus mirabilis</i> ATCC® 7002™
<b>MSP018</b>	MicroSwab <i>Proteus hauseri</i> ATCC® 13315™
<b>MSP019</b>	MicroSwab <i>Proteus vulgaris</i> ATCC® 49132™
<b>MSP020</b>	MicroSwab <i>Proteus vulgaris</i> ATCC® 6380™
<b>MSP021</b>	MicroSwab <i>Proteus vulgaris</i> ATCC® 8427™
<b>MSP022</b>	MicroSwab <i>Providencia alcalifaciens</i> ATCC® 51902™
<b>MSP023</b>	MicroSwab <i>Providencia stuartii</i> ATCC® 33672™
<b>MSP024</b>	MicroSwab <i>Providencia stuartii</i> ATCC® 49809™
<b>MSP025</b>	MicroSwab <i>Prevotella loeschii</i> ATCC® 15930™
<b>MSP026</b>	MicroSwab <i>Pseudomonas aeruginosa</i> WDCM 00024-ATCC® 10145™
<b>MSP027</b>	MicroSwab <i>Pseudomonas aeruginosa</i> ATCC® 25619™
<b>MSP028</b>	MicroSwab <i>Pseudomonas aeruginosa</i> WDCM 00025-ATCC® 27853™
<b>MSP029</b>	MicroSwab <i>Pseudomonas aeruginosa</i> ATCC® 35422™
<b>MSP030</b>	MicroSwab <i>Pseudomonas aeruginosa</i> WDCM 00026-ATCC® 9027™
<b>MSP032</b>	MicroSwab <i>Pseudomonas fluorescens</i> WDCM 00115-ATCC® 13525™
<b>MSP033</b>	MicroSwab <i>Pseudomonas putida</i> ATCC® 49128™
<b>MSP034</b>	MicroSwab <i>Pseudomonas stutzeri</i> ATCC® 17588™
<b>MSP035</b>	MicroSwab <i>Proteus mirabilis</i> ATCC® 25933™
<b>MSP036</b>	MicroSwab <i>Pseudomonas aeruginosa</i> ATCC® 19429™
<b>MSP038</b>	MicroSwab <i>Pseudomonas aeruginosa</i> ATCC® 15442™
<b>MSP040</b>	MicroSwab <i>Parabacteroides distasonis</i> ATCC® BAA1295™
<b>MSP041</b>	MicroSwab <i>Phytophthora sojae</i> ATCC® 48069™
<b>MSP043</b>	MicroSwab <i>Prevotella melaninogenica</i> ATCC® 25845™
<b>MSP044</b>	MicroSwab <i>Pseudomonas aeruginosa</i> WDCM 00114-ATCC® 25668™
<b>MSP045</b>	MicroSwab <i>Pseudomonas</i> sp. ATCC® 19144™
<b>MSR001</b>	MicroSwab <i>Rhizopus stolonifer</i> ATCC® 14037™
<b>MSR002</b>	MicroSwab <i>Rhodococcus equi</i> WDCM 00028- ATCC® 6939™
<b>MSS001</b>	MicroSwab <i>Saccharomyces cerevisiae</i> ATCC® 32701™
<b>MSS002</b>	MicroSwab <i>Saccharomyces cerevisiae</i> ATCC® 4098™



© The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalogue marks are trademarks of ATCC. Meccanti S.A.R.L. Sp. z o. o. is licensed to use these trademarks and sell products derived from ATCC cultures. Look for the ATCC Licensed Derivative Emblem for products derived from ATCC cultures.

# MICROSWAB

code	product name
<b>MSS003</b>	MicroSwab <i>Salmonella abaetetuba</i> ATCC® 35640™
<b>MSS004</b>	MicroSwab <i>Salmonella anatum</i> ATCC® 9270™
<b>MSS005</b>	MicroSwab <i>Salmonella arizonae</i> ATCC® 13314™
<b>MSS006</b>	MicroSwab <i>Salmonella enteritidis</i> WDCM 00030-ATCC® 13076™
<b>MSS008</b>	MicroSwab <i>Salmonella diarizonae</i> ATCC® 29226™
<b>MSS009</b>	MicroSwab <i>Salmonella typhimurium</i> ATCC® 29631™
<b>MSS011</b>	MicroSwab <i>Salmonella schottmuelleri</i> ATCC® 8759™
<b>MSS014</b>	MicroSwab <i>Salmonella kunzendorf</i> ATCC® 12011™
<b>MSS015</b>	MicroSwab <i>Salmonella paratyphi</i> ATCC® 9150™
<b>MSS016</b>	MicroSwab <i>Salmonella tallahassee</i> ATCC® 12002™
<b>MSS018</b>	MicroSwab <i>Salmonella typhimurium</i> WDCM 00031-ATCC® 14028™
<b>MSS019</b>	MicroSwab <i>Salmonella vellore</i> ATCC® 15611™
<b>MSS020</b>	MicroSwab <i>Serratia liquefaciens</i> ATCC® 27592™
<b>MSS021</b>	MicroSwab <i>Serratia marcescens</i> ATCC® 13880™
<b>MSS022</b>	MicroSwab <i>Serratia marcescens</i> ATCC® 14756™
<b>MSS023</b>	MicroSwab <i>Serratia marcescens</i> ATCC® 8100™
<b>MSS025</b>	MicroSwab <i>Serratia odorifera</i> ATCC® 33077™
<b>MSS026</b>	MicroSwab <i>Serratia rubidaea</i> ATCC® 33670™
<b>MSS027</b>	MicroSwab <i>Shewanella putrefaciens</i> ATCC® 49138™
<b>MSS028</b>	MicroSwab <i>Shewanella putrefaciens</i> ATCC® 8071™
<b>MSS029</b>	MicroSwab <i>Shigella boydii</i> ATCC® 9207™
<b>MSS031</b>	MicroSwab <i>Shigella flexneri</i> WDCM 00126-ATCC® 12022™
<b>MSS032</b>	MicroSwab <i>Shigella sonnei</i> ATCC® 25931™
<b>MSS033</b>	MicroSwab <i>Shigella sonnei</i> ATCC® 9290™
<b>MSS034</b>	MicroSwab <i>Sphingobacterium multivorum</i> ATCC® 35656™
<b>MSS035</b>	MicroSwab <i>Sphingomonas trueperi</i> ATCC® 12417™
<b>MSS036</b>	MicroSwab <i>Staphylococcus aureus</i> ATCC® 12598™
<b>MSS037</b>	MicroSwab <i>Staphylococcus aureus</i> ATCC® 12600™
<b>MSS038</b>	MicroSwab <i>Staphylococcus aureus</i> ATCC® 25904™
<b>MSS039</b>	MicroSwab <i>Staphylococcus aureus</i> WDCM 00034-ATCC® 25923™
<b>MSS040</b>	MicroSwab <i>Staphylococcus aureus</i> WDCM 00131-ATCC® 29213™
<b>MSS041</b>	MicroSwab <i>Staphylococcus aureus</i> ATCC® 31153™
<b>MSS042</b>	MicroSwab <i>Staphylococcus aureus</i> ATCC® 33591™



© The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalogue marks are trademarks of ATCC. Meccanti S.p.A., Sp. z o. o. is licensed to use these trademarks and sell products derived from ATCC cultures. Look for the ATCC Licensed Derivative Emblem for products derived from ATCC cultures.

# MICROSWAB

code	product name
<b>MSS043</b>	MicroSwab Staphylococcus aureus ATCC® 33862™
<b>MSS044</b>	MicroSwab Staphylococcus aureus WDCM 00211-ATCC® 43300™
<b>MSS045</b>	MicroSwab Staphylococcus aureus ATCC® 49476™
<b>MSS046</b>	MicroSwab Staphylococcus aureus ATCC® 51153™
<b>MSS047</b>	MicroSwab Staphylococcus aureus WDCM 00032,WDCM 00193-ATCC® 6538™
<b>MSS048</b>	MicroSwab Staphylococcus aureus WDCM 00033,WDCM 00193,WDCM 00195-ATCC® 6538P™
<b>MSS049</b>	MicroSwab Staphylococcus capitis ATCC® 35661™
<b>MSS050</b>	MicroSwab Macrococcus caseolyticus ATCC® 35662™
<b>MSS051</b>	MicroSwab Staphylococcus epidermidis WDCM 00036-ATCC® 12228™
<b>MSS052</b>	MicroSwab Staphylococcus epidermidis WDCM 00132-ATCC® 14990™
<b>MSS053</b>	MicroSwab Staphylococcus epidermidis ATCC® 29887™
<b>MSS054</b>	MicroSwab Staphylococcus epidermidis ATCC® 49134™
<b>MSS055</b>	MicroSwab Staphylococcus hominis ATCC® 27844™
<b>MSS056</b>	MicroSwab Staphylococcus lentus ATCC® 700403™
<b>MSS057</b>	MicroSwab Staphylococcus saprophyticus WDCM 00159-ATCC® 15305™
<b>MSS058</b>	MicroSwab Staphylococcus saprophyticus ATCC® 35552™
<b>MSS059</b>	MicroSwab Staphylococcus saprophyticus ATCC® 49453™
<b>MSS060</b>	MicroSwab Staphylococcus saprophyticus ATCC® 49907™
<b>MSS061</b>	MicroSwab Staphylococcus sciuri ATCC® 29060™
<b>MSS062</b>	MicroSwab Staphylococcus xylosus ATCC® 29967™
<b>MSS063</b>	MicroSwab Staphylococcus xylosus ATCC® 29971™
<b>MSS064</b>	MicroSwab Staphylococcus xylosus ATCC® 35663™
<b>MSS065</b>	MicroSwab Staphylococcus xylosus ATCC® 49148™
<b>MSS066</b>	MicroSwab Staphylococcus xylosus ATCC® 700404™
<b>MSS067</b>	MicroSwab Stenotrophomonas maltophilia ATCC® 13637™
<b>MSS068</b>	MicroSwab Stenotrophomonas maltophilia ATCC® 51331™
<b>MSS069</b>	MicroSwab Streptococcus agalactiae ATCC® 13813™
<b>MSS070</b>	MicroSwab Streptococcus agalactiae ATCC® 27956™
<b>MSS072</b>	MicroSwab Streptococcus anginosus ATCC® 33397™
<b>MSS073</b>	MicroSwab Streptococcus dysgalactiae ATCC® 35666™
<b>MSS074</b>	MicroSwab Streptococcus pasteurianus ATCC® 49133™



© The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalogue marks are trademarks of ATCC. Mecconti S.A.R.L. Sp. z o. o. is licensed to use these trademarks and sell products derived from ATCC cultures. Look for the ATCC Licensed Derivative Emblem for products derived from ATCC cultures.

# MICROSWAB

code	product name
<b>MSS075</b>	MicroSwab Streptococcus gallolyticus ATCC® 49475™
<b>MSS076</b>	MicroSwab Streptococcus bovis ATCC® 35034™
<b>MSS077</b>	MicroSwab Streptococcus gallolyticus ATCC® 49147™
<b>MSS078</b>	MicroSwab Streptococcus gallolyticus ATCC® 9809™
<b>MSS079</b>	MicroSwab Streptococcus equi subsp.equi ATCC® 9528™
<b>MSS080</b>	MicroSwab Streptococcus oralis ATCC® 6249™
<b>MSS081</b>	MicroSwab Streptococcus mutans ATCC® 35668™
<b>MSS082</b>	MicroSwab Streptococcus oralis ATCC® 9811™
<b>MSS083</b>	MicroSwab Streptococcus pneumoniae ATCC® 27336™
<b>MSS084</b>	MicroSwab Streptococcus pneumoniae ATCC® 49136™
<b>MSS085</b>	MicroSwab Streptococcus pneumoniae ATCC® 49150™
<b>MSS086</b>	MicroSwab Streptococcus pneumoniae ATCC® 49619™
<b>MSS087</b>	MicroSwab Streptococcus pneumoniae ATCC® 6303™
<b>MSS088</b>	MicroSwab Streptococcus pneumoniae ATCC® 6305™
<b>MSS089</b>	MicroSwab Streptococcus pyogenes ATCC® 12384™
<b>MSS090</b>	MicroSwab Streptococcus pyogenes ATCC® 19615™
<b>MSS091</b>	MicroSwab Streptococcus sanguis ATCC® 10556™
<b>MSS092</b>	MicroSwab Streptococcus agalactiae ATCC® 12386™
<b>MSS093</b>	MicroSwab Streptococcus dysgalactiae ATCC® 12388™
<b>MSS094</b>	MicroSwab Streptococcus sp. (Group F) ATCC® 12392™
<b>MSS095</b>	MicroSwab Streptococcus dysgalactiae ATCC® 12394™
<b>MSS096</b>	MicroSwab Streptococcus sp. (Group D) ATCC® 27284™
<b>MSS097</b>	MicroSwab Streptococcus uberis ATCC® 700407™
<b>MSS099</b>	MicroSwab Streptomyces somialensis ATCC® 33201™
<b>MSS100</b>	MicroSwab Salmonella abony WDCM 00029-NCTC® 6017™
<b>MSS101</b>	MicroSwab Saccharomyces cerevisiae WDCM 00058-ATCC® 9763™
<b>MSS102</b>	MicroSwab Salmonella choleraesuis ATCC® 10708™
<b>MSS103</b>	MicroSwab Salmonella typhimurium WDCM 00121-ATCC® 13311™
<b>MSS104</b>	MicroSwab Salmonella pullorum ATCC® 19945™
<b>MSS105</b>	MicroSwab Staphylococcus aureus WDCM 00035-ATCC® 9144™
<b>MSS106</b>	MicroSwab Staphylococcus aureus ATCC® BAA 1026™
<b>MSS107</b>	MicroSwab Staphylococcus aureus ATCC® BAA 977™
<b>MSS108</b>	MicroSwab Staphylococcus aureus ATCC® BAA 976™
<b>MSS112</b>	MicroSwab Stenotrophomonas maltophilia ATCC® 17666™



The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalogue marks are trademarks of ATCC. Meccanti S.p.A. is licensed to use these trademarks and sell products derived from ATCC cultures. Look for the ATCC Licensed Derivative Emblem for products derived from ATCC cultures.

# MICROSWAB

code	product name
<b>MSS113</b>	MicroSwab Salmonella Group D, not typhi ATCC® CL BF-SD™
<b>MSS114</b>	MicroSwab Salmonella sp. ATCC® 35664™
<b>MSS115</b>	MicroSwab Serratia marcescens ATCC® 27857™
<b>MSS116</b>	MicroSwab Serratia marcescens ATCC® CLD1™
<b>MSS118</b>	MicroSwab Staphylococcus intermedius ATCC® 29663™
<b>MSS119</b>	MicroSwab Staphylococcus pseudintermedius ATCC® 49444™
<b>MSS120</b>	MicroSwab Staphylococcus saprophyticus ATCC® BAA750™
<b>MSS121</b>	MicroSwab Staphylococcus sciuri ATCC® 29061™
<b>MSS122</b>	MicroSwab Streptococcus agalactiae ATCC® CL 810™
<b>MSS123</b>	MicroSwab Streptococcus bovis WDCM 00133-ATCC® 33317™
<b>MSS124</b>	MicroSwab Streptococcus equi subsp. zooepidemicus ATCC® 700400™
<b>MSS125</b>	MicroSwab Streptococcus mutans ATCC® 25175™
<b>MSS126</b>	MicroSwab Streptococcus salivarius ATCC® 13419™
<b>MSS127</b>	MicroSwab Salmonella enterica ATCC® 12325™
<b>MST001</b>	MicroSwab Trichophyton mentagrophytes ATCC® 9533™
<b>MST002</b>	MicroSwab Trichophyton rubrum ATCC® 28188™
<b>MST003</b>	MicroSwab Tannerella forsythensis ATCC® 700191™
<b>MSV001</b>	MicroSwab Vibrio parahaemolyticus WDCM 00037-ATCC® 17802™
<b>MSY001</b>	MicroSwab Yarrowia lipolytica ATCC® 9773™
<b>MSY002</b>	MicroSwab Yersinia enterocolitica WDCM 00160-ATCC® 23715™
<b>MSY003</b>	MicroSwab Yersinia enterocolitica WDCM 00038-ATCC® 9610™
<b>MSY004</b>	MicroSwab Yersinia kristensenii ATCC® 33639™
<b>MSZ001</b>	MicroSwab Zygosaccharomyces rouxii ATCC® 28253™

# MICROSWABPLUS

code	product name
<b>MS2A001</b>	MicroSwabPlus Acinetobacter baumannii ATCC® 19606™
<b>MS2A009</b>	MicroSwabPlus Aerococcus viridans ATCC® 700406™
<b>MS2A012</b>	MicroSwabPlus Aeromonas hydrophila WDCM 00063-ATCC® 7966™
<b>MS2A017</b>	MicroSwabPlus Aspergillus brasiliensis WDCM 00053-ATCC® 16404™
<b>MS2B001</b>	MicroSwabPlus Bacillus cereus ATCC® 10876™



The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalogue marks are trademarks of ATCC. Mecconti S.p.A.R.L. is licensed to use these trademarks and sell products derived from ATCC cultures. Look for the ATCC Licensed Derivative Emblem for products derived from ATCC cultures.

## MICROORGANISMS

LYOPHILIZED REFERENCE ORGANISMS

MECCONTI S.p.A.R.L. SP.Z O.O.

# MICROSWABPLUS

code	product name
<b>MS2B002</b>	MicroSwabPlus <i>Bacillus cereus</i> WDCM 00001-ATCC® 11778™
<b>MS2B006</b>	MicroSwabPlus <i>Geobacillus stearothermophilus</i> ATCC® 7953™
<b>MS2B007</b>	MicroSwabPlus <i>Bacillus subtilis</i> subsp. <i>spizizenii</i> WDCM 00003-ATCC® 6633™
<b>MS2B010</b>	MicroSwabPlus <i>Bacteroides fragilis</i> ATCC® 25285™
<b>MS2B014</b>	MicroSwabPlus <i>Bacteroides ureolyticus</i> ATCC® 33387™
<b>MS2B016</b>	MicroSwabPlus <i>Bacteroides vulgatus</i> ATCC® 8482™
<b>MS2B018</b>	MicroSwabPlus <i>Bordetella bronchiseptica</i> ATCC® 10580™
<b>MS2B019</b>	MicroSwabPlus <i>Bordetella bronchiseptica</i> ATCC® 4617™
<b>MS2B022</b>	MicroSwabPlus <i>Moraxella catarrhalis</i> ATCC® 25238™
<b>MS2B026</b>	MicroSwabPlus <i>Brevundimonas diminuta</i> ATCC® 19146™
<b>MS2B027</b>	MicroSwabPlus <i>Burkholderia cepacia</i> ATCC® 25416™
<b>MS2B028</b>	MicroSwabPlus <i>Burkholderia cepacia</i> ATCC® 25608™
<b>MS2B031</b>	MicroSwabPlus <i>Bacillus pumilus</i> ATCC® 14884™
<b>MS2C002</b>	MicroSwabPlus <i>Campylobacter jejuni</i> WDCM 00156-ATCC® 29428™
<b>MS2C003</b>	MicroSwabPlus <i>Campylobacter jejuni</i> WDCM 00005-ATCC® 33291™
<b>MS2C005</b>	MicroSwabPlus <i>Candida albicans</i> WDCM 00054-ATCC® 10231™
<b>MS2C009</b>	MicroSwabPlus <i>Candida glabrata</i> ATCC® 15126™
<b>MS2C022</b>	MicroSwabPlus <i>Citrobacter freundii</i> ATCC® 8090™
<b>MS2C023</b>	MicroSwabPlus <i>Clostridium difficile</i> ATCC® 9689™
<b>MS2C026</b>	MicroSwabPlus <i>Clostridium perfringens</i> WDCM 00007-ATCC® 13124™
<b>MS2C029</b>	MicroSwabPlus <i>Clostridium sporogenes</i> ATCC® 11437™
<b>MS2C030</b>	MicroSwabPlus <i>Clostridium sporogenes</i> WDCM 00008-ATCC® 19404™
<b>MS2C049</b>	MicroSwabPlus <i>Candida albicans</i> ATCC® 90028™
<b>MS2C050</b>	MicroSwabPlus <i>Campylobacter coli</i> WDCM 00004-ATCC® 43478™
<b>MS2C051</b>	MicroSwabPlus <i>Citrobacter freundii</i> WDCM 00006-ATCC® 43864™
<b>MS2E001</b>	MicroSwabPlus <i>Enterobacter aerogenes</i> WDCM 00175-ATCC® 13048™
<b>MS2E003</b>	MicroSwabPlus <i>Enterobacter cloacae</i> WDCM 00083-ATCC® 13047™
<b>MS2E004</b>	MicroSwabPlus <i>Enterobacter cloacae</i> WDCM 00082-ATCC® 23355™
<b>MS2E012</b>	MicroSwabPlus <i>Enterococcus durans</i> ATCC® 49479™
<b>MS2E013</b>	MicroSwabPlus <i>Enterococcus durans</i> ATCC® 6056™



© The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalogue marks are trademarks of ATCC. Meccanti S.p.A. is licensed to use these trademarks and sell products derived from ATCC cultures. Look for the ATCC Licensed Derivative Emblem for products derived from ATCC cultures.

# MICROSWABPLUS

code	product name
<b>MS2E014</b>	MicroSwabPlus Enterococcus faecalis WDCM 00009-ATCC® 19433™
<b>MS2E015</b>	MicroSwabPlus Enterococcus faecalis WDCM 00087-ATCC® 29212™
<b>MS2E016</b>	MicroSwabPlus Enterococcus faecalis ATCC® 49149™
<b>MS2E017</b>	MicroSwabPlus Enterococcus faecalis ATCC® 51299™
<b>MS2E018</b>	MicroSwabPlus Enterococcus faecium ATCC® 35667™
<b>MS2E022</b>	MicroSwabPlus Escherichia coli ATCC® 11303™
<b>MS2E023</b>	MicroSwabPlus Escherichia coli WDCM 00090-ATCC® 11775™
<b>MS2E025</b>	MicroSwabPlus Escherichia coli WDCM 00013-ATCC® 25922™
<b>MS2E027</b>	MicroSwabPlus Escherichia coli ATCC® 35218™
<b>MS2E032</b>	MicroSwabPlus Escherichia coli WDCM 00012,WDCM 00196-ATCC® 8739™
<b>MS2E034</b>	MicroSwabPlus Cronobacter muytjensii WDCM 00213- ATCC® 51329™
<b>MS2E036</b>	MicroSwabPlus Enterococcus hirae WDCM 00011-ATCC® 10541™
<b>MS2E037</b>	MicroSwabPlus Escherichia coli ATCC® 10536™
<b>MS2E039</b>	MicroSwabPlus Enterococcus casseliflavus ATCC® 700327™
<b>MS2E040</b>	MicroSwabPlus Enterococcus faecium WDCM 00177-ATCC® 6057™
<b>MS2E048</b>	MicroSwabPlus Enterobacter hormaechei ATCC® 700323™
<b>MS2F002</b>	MicroSwabPlus Fusobacterium nucleatum ATCC® 25586™
<b>MS2H003</b>	MicroSwabPlus Haemophilus influenzae ATCC® 10211™
<b>MS2H009</b>	MicroSwabPlus Haemophilus influenzae ATCC® 49247™
<b>MS2H010</b>	MicroSwabPlus Haemophilus influenzae ATCC® 49766™
<b>MS2H011</b>	MicroSwabPlus Haemophilus influenzae ATCC® 9006™
<b>MS2H015</b>	MicroSwabPlus Haemophilus parainfluenzae ATCC® 7901™
<b>MS2K001</b>	MicroSwabPlus Klebsiella oxytoca ATCC® 43086™
<b>MS2K002</b>	MicroSwabPlus Klebsiella oxytoca ATCC® 49131™
<b>MS2K003</b>	MicroSwabPlus Klebsiella pneumoniae ATCC® 10031™
<b>MS2K004</b>	MicroSwabPlus Klebsiella pneumoniae WDCM 00097-ATCC® 13883™
<b>MS2K006</b>	MicroSwabPlus Klebsiella pneumoniae ATCC® 27736™
<b>MS2K008</b>	MicroSwabPlus Klebsiella pneumoniae ATCC® 35657™
<b>MS2K010</b>	MicroSwabPlus Klebsiella pneumoniae ATCC® 9997™
<b>MS2K012</b>	MicroSwabPlus Klebsiella oxytoca ATCC® 700324™
<b>MS2K013</b>	MicroSwabPlus Klebsiella pneumoniae ATCC® BAA1705™



© The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalogue marks are trademarks of ATCC. Mecconti S.A.R.L. Sp. z o. o. is licensed to use these trademarks and sell products derived from ATCC cultures. Look for the ATCC Licensed Derivative Emblem for products derived from ATCC cultures.

# MICROSWABPLUS

code	product name
<b>MS2K014</b>	MicroSwabPlus Klebsiella pneumoniae ATCC® BAA1706™
<b>MS2K016</b>	MicroSwabPlus Klebsiella variicola WDCM 00206-ATCC® 31488™
<b>MS2L001</b>	MicroSwabPlus Lactobacillus acidophilus WDCM 00098-ATCC® 4356™
<b>MS2L003</b>	MicroSwabPlus Lactobacillus rhamnosus WDCM 00101-ATCC® 7469™
<b>MS2L005</b>	MicroSwabPlus Leclercia adecarboxylata ATCC® 23216™
<b>MS2L006</b>	MicroSwabPlus Legionella pneumophila WDCM 00107-ATCC® 33152™
<b>MS2L008</b>	MicroSwabPlus Listeria innocua WDCM 00017-ATCC® 33090™
<b>MS2L010</b>	MicroSwabPlus Listeria ivanovii WDCM 00018-ATCC® 19119™
<b>MS2L011</b>	MicroSwabPlus Listeria monocytogenes WDCM 00020-ATCC® 19111™
<b>MS2L012</b>	MicroSwabPlus Listeria monocytogenes ATCC® 7644™
<b>MS2L013</b>	MicroSwabPlus Listeria monocytogenes ATCC® 19115™
<b>MS2L023</b>	MicroSwabPlus Listeria monocytogenes WDCM 00021-ATCC® 13932™
<b>MS2M001</b>	MicroSwabPlus Micrococcus luteus WDCM 00111-ATCC® 4698™
<b>MS2M004</b>	MicroSwabPlus Kocuria rhizophila ATCC® 9341™
<b>MS2M005</b>	MicroSwabPlus Micrococcus sp. ATCC® 700405™
<b>MS2M009</b>	MicroSwabPlus Micrococcus luteus ATCC® 10240™
<b>MS2N003</b>	MicroSwabPlus Neisseria gonorrhoeae ATCC® 43069™
<b>MS2N005</b>	MicroSwabPlus Neisseria gonorrhoeae ATCC® 49226™
<b>MS2N006</b>	MicroSwabPlus Neisseria gonorrhoeae ATCC® 49981™
<b>MS2N007</b>	MicroSwabPlus Neisseria lactamica ATCC® 23970™
<b>MS2N012</b>	MicroSwabPlus Neisseria meningitidis ATCC® 13102™(serogrup C)
<b>MS2P002</b>	MicroSwabPlus Penicillium chrysogenum ATCC® 10106™
<b>MS2P003</b>	MicroSwabPlus Peptostreptococcus anaerobius ATCC® 27337™
<b>MS2P006</b>	MicroSwabPlus Micromonas micros ATCC® 33270™
<b>MS2P008</b>	MicroSwabPlus Plesiomonas shigelloides ATCC® 51903™
<b>MS2P011</b>	MicroSwabPlus Propionibacterium acnes ATCC® 11827™



© The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalogue marks are trademarks of ATCC. Meccanti S.p.A., Sp. z o. o. is licensed to use these trademarks and sell products derived from ATCC cultures. Look for the ATCC Licensed Derivative Emblem for products derived from ATCC cultures.

# MICROSWABPLUS

code	product name
<b>MS2P012</b>	MicroSwabPlus <i>Proteus mirabilis</i> ATCC® 12453™
<b>MS2P013</b>	MicroSwabPlus <i>Proteus mirabilis</i> ATCC® 29245™
<b>MS2P014</b>	MicroSwabPlus <i>Proteus mirabilis</i> WDCM 00023-ATCC® 29906™
<b>MS2P015</b>	MicroSwabPlus <i>Proteus mirabilis</i> ATCC® 35659™
<b>MS2P017</b>	MicroSwabPlus <i>Proteus mirabilis</i> ATCC® 7002™
<b>MS2P018</b>	MicroSwabPlus <i>Proteus hauseri</i> ATCC® 13315™
<b>MS2P020</b>	MicroSwabPlus <i>Proteus vulgaris</i> ATCC® 6380™
<b>MS2P021</b>	MicroSwabPlus <i>Proteus vulgaris</i> ATCC® 8427™
<b>MS2P022</b>	MicroSwabPlus <i>Providencia alcalifaciens</i> ATCC® 51902™
<b>MS2P026</b>	MicroSwabPlus <i>Pseudomonas aeruginosa</i> WDCM 00024-ATCC® 10145™
<b>MS2P028</b>	MicroSwabPlus <i>Pseudomonas aeruginosa</i> WDCM 00025-ATCC® 27853™
<b>MS2P030</b>	MicroSwabPlus <i>Pseudomonas aeruginosa</i> WDCM 00026-ATCC® 9027™
<b>MS2P032</b>	MicroSwabPlus <i>Pseudomonas fluorescens</i> WDCM 00115-ATCC® 13525™
<b>MS2P035</b>	MicroSwabPlus <i>Proteus mirabilis</i> ATCC® 25933™
<b>MS2P038</b>	MicroSwabPlus <i>Pseudomonas aeruginosa</i> ATCC® 15442™
<b>MS2R002</b>	MicroSwabPlus <i>Rhodococcus equi</i> WDCM 00028- ATCC® 6939™
<b>MS2S001</b>	MicroSwabPlus <i>Saccharomyces cerevisiae</i> ATCC® 32701™
<b>MS2S004</b>	MicroSwabPlus <i>Salmonella anatum</i> ATCC® 9270™
<b>MS2S005</b>	MicroSwabPlus <i>Salmonella arizonaee</i> ATCC® 13314™
<b>MS2S006</b>	MicroSwabPlus <i>Salmonella enteritidis</i> WDCM 00030-ATCC® 13076™
<b>MS2S011</b>	MicroSwabPlus <i>Salmonella schottmuelleri</i> ATCC® 8759™
<b>MS2S015</b>	MicroSwabPlus <i>Salmonella paratyphi</i> ATCC® 9150™
<b>MS2S018</b>	MicroSwabPlus <i>Salmonella typhimurium</i> WDCM 00031-ATCC® 14028™
<b>MS2S020</b>	MicroSwabPlus <i>Serratia liquefaciens</i> ATCC® 27592™
<b>MS2S021</b>	MicroSwabPlus <i>Serratia marcescens</i> ATCC® 13880™
<b>MS2S023</b>	MicroSwabPlus <i>Serratia marcescens</i> ATCC® 8100™
<b>MS2S029</b>	MicroSwabPlus <i>Shigella boydii</i> ATCC® 9207™
<b>MS2S031</b>	MicroSwabPlus <i>Shigella flexneri</i> WDCM 00126-ATCC® 12022™
<b>MS2S032</b>	MicroSwabPlus <i>Shigella sonnei</i> ATCC® 25931™



© The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalogue marks are trademarks of ATCC. Mecconti S.p.A.R.L. Sp. z o. o. is licensed to use these trademarks and sell products derived from ATCC cultures. Look for the ATCC Licensed Derivative Emblem for products derived from ATCC cultures.

# MICROSWABPLUS

code	product name
<b>MS2S033</b>	MicroSwabPlus Shigella sonnei ATCC® 9290™
<b>MS2S039</b>	MicroSwabPlus Staphylococcus aureus WDCM 00034-ATCC® 25923™
<b>MS2S040</b>	MicroSwabPlus Staphylococcus aureus WDCM 00131-ATCC® 29213™
<b>MS2S042</b>	MicroSwabPlus Staphylococcus aureus ATCC® 33591™
<b>MS2S044</b>	MicroSwabPlus Staphylococcus aureus WDCM 00211-ATCC® 43300™
<b>MS2S047</b>	MicroSwabPlus Staphylococcus aureus WDCM 00032,WDCM 00193-ATCC® 6538™
<b>MS2S048</b>	MicroSwabPlus Staphylococcus aureus WDCM 00033,WDCM 00193,WDCM 00195-ATCC® 6538P™
<b>MS2S049</b>	MicroSwabPlus Staphylococcus capitis ATCC® 35661™
<b>MS2S051</b>	MicroSwabPlus Staphylococcus epidermidis WDCM 00036-ATTC® 12228™
<b>MS2S052</b>	MicroSwabPlus Staphylococcus epidermidis WDCM 00132-ATCC® 14990™
<b>MS2S054</b>	MicroSwabPlus Staphylococcus epidermidis ATCC® 49134™
<b>MS2S057</b>	MicroSwabPlus Staphylococcus saprophyticus WDCM 00159-ATCC® 15305™
<b>MS2S059</b>	MicroSwabPlus Staphylococcus saprophyticus ATCC® 49453™
<b>MS2S066</b>	MicroSwabPlus Staphylococcus xylosus ATCC® 700404™
<b>MS2S069</b>	MicroSwabPlus Streptococcus agalactiae ATCC® 13813™
<b>MS2S072</b>	MicroSwabPlus Streptococcus anginosus ATCC® 33397™
<b>MS2S081</b>	MicroSwabPlus Streptococcus mutans ATCC® 35668™
<b>MS2S086</b>	MicroSwabPlus Streptococcus pneumoniae ATCC® 49619™
<b>MS2S087</b>	MicroSwabPlus Streptococcus pneumoniae ATCC® 6303™
<b>MS2S090</b>	MicroSwabPlus Streptococcus pyogenes ATCC® 19615™
<b>MS2S097</b>	MicroSwabPlus Streptococcus uberis ATCC® 700407™
<b>MS2S100</b>	MicroSwabPlus Salmonella abony WDCM 00029-NCTC® 6017™
<b>MS2S101</b>	MicroSwabPlus Saccharomyces cerevisiae WDCM 00058-ATCC® 9763™
<b>MS2S102</b>	MicroSwabPlus Salmonella choleraesuis ATCC® 10708™
<b>MS2S103</b>	MicroSwabPlus Salmonella typhimurium WDCM 00121-ATCC® 13311™
<b>MS2S108</b>	MicroSwabPlus Staphylococcus aureus ATCC® BAA 976™
<b>MS2S112</b>	MicroSwabPlus Stenotrophomonas maltophilia ATCC® 17666™



The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalogue marks are trademarks of ATCC. Meccanti S.p.A. is licensed to use these trademarks and sell products derived from ATCC cultures. Look for the ATCC Licensed Derivative Emblem for products derived from ATCC cultures.

# MICROSWABPLUS

code	product name
<b>MS2S124</b>	MicroSwabPlus Streptococcus equi subsp. zooepidemicus ATCC® 700400™
<b>MS2T001</b>	MicroSwabPlus Trichophyton mentagrophytes ATCC® 9533™
<b>MS2V001</b>	MicroSwabPlus Vibrio parahaemolyticus WDCM 00037-ATCC® 17802™
<b>MS2Y002</b>	MicroSwabPlus Yersinia enterocolitica WDCM 00160-ATCC® 23715™
<b>MS2Y003</b>	MicroSwabPlus Yersinia enterocolitica WDCM 00038-ATCC® 9610™
<b>MS2Y004</b>	MicroSwabPlus Yersinia kristensenii ATCC® 33639™
<b>MS2Z001</b>	MicroSwabPlus Zygosaccharomyces rouxii ATCC® 28253™

# MICROPELLETRT

code	product name
<b>MPRTA017</b>	MicroPelletRT Aspergillus brasiliensis WDCM 00053-ATCC® 16404™
<b>MPRTB001</b>	MicroPelletRT Bacillus cereus ATCC® 10876™
<b>MPRTB006</b>	MicroPelletRT Geobacillus stearothermophilus ATCC® 7953™
<b>MPRTB007</b>	MicroPelletRT Bacillus subtilis subsp. Spizizenii WDCM 00003-ATCC® 6633™
<b>MPRTB026</b>	MicroPelletRT Brevundimonas diminuta ATCC® 19146
<b>MPRTB027</b>	MicroPelletRT Burkholderia cepacia ATCC® 25416
<b>MPRTC002</b>	MicroPelletRT Candida albicans WDCM 00055-ATCC® 2091™
<b>MPRTC005</b>	MicroPelletRT Candida albicans WDCM 00054-ATCC® 10231™
<b>MPRTC026</b>	MicroPelletRT Clostridium perfringens WDCM 00007-ATCC® 13124™
<b>MPRTC029</b>	MicroPelletRT Clostridium sporogenes ATCC® 11437™
<b>MPRTC030</b>	MicroPelletRT Clostridium sporogenes WDCM 00008-ATCC® 19404™
<b>MPRTE001</b>	MicroPelletRT Klebsiella aerogenes WDCM 00175-ATCC® 13048™
<b>MPRTE015</b>	MicroPelletRT Enterococcus faecalis WDCM 00087-ATCC® 29212™
<b>MPRTE032</b>	MicroPelletRT Escherichia coli WDCM 00012, WDCM 00196-ATCC® 8739™



The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalogue marks are trademarks of ATCC. Meccanti S.p.A.R.L. Sp. z o. o. is licensed to use these trademarks and sell products derived from ATCC cultures. Look for the ATCC Licensed Derivative Emblem for products derived from ATCC cultures.

## MICROORGANISMS

LYOPHILIZED REFERENCE ORGANISMS

MECCANTI S.p.A.R.L. Sp. z o.o.

# MICROPELLET<sup>RT</sup>

code	product name
<b>MPRTK013</b>	MicroPelletRT <i>Kocuria rhizophila</i> ATCC® 9341™
<b>MPRTL013</b>	MicroPelletRT <i>Listeria monocytogenes</i> ATCC® 19115™
<b>MPRTL016</b>	MicroPelletRT <i>Lactobacillus fermentum</i> ATCC® 9338™
<b>MPRTP028</b>	MicroPelletRT <i>Pseudomonas aeruginosa</i> WDCM 00025-ATCC® 27853™
<b>MPRTP030</b>	MicroPelletRT <i>Pseudomonas aeruginosa</i> WDCM 00026-ATCC® 9027™
<b>MPRTS003</b>	MicroPelletRT <i>Salmonella enterica</i> subsp. <i>Enterica</i> serovar <i>Abaetetuba</i> ATCC® 35640™
<b>MPRTS018</b>	MicroPelletRT <i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Typhimurium</i> WDCM 00031-ATCC® 14028™
<b>MPRTS039</b>	MicroPelletRT <i>Staphylococcus aureus</i> subsp. <i>aureus</i> WDCM 00034-ATCC® 25923™
<b>MPRTS047</b>	MicroPelletRT <i>Staphylococcus aureus</i> subsp. <i>aureus</i> WDCM 00032,WDCM 00193 -ATCC® 6538™
<b>MPRTS048</b>	MicroPelletRT <i>Staphylococcus aureus</i> subsp. <i>aureus</i> WDCM 00033,WDCM 00193,WDCM 00195-ATCC® 6538P™
<b>MPRTS051</b>	MicroPelletRT <i>Staphylococcus epidermidis</i> WDCM 00036-ATCC® 12228™
<b>MPRTS090</b>	MicroPelletRT <i>Streptococcus pyogenes</i> ATCC® 19615™
<b>MPRTS100</b>	MicroPelletRT <i>Salmonella abony</i> WDCM 00029-NCTC® 6017™
<b>MPRTS102</b>	MicroPelletRT <i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Choleraesuis</i> ATCC® 10708™
<b>MPRTS103</b>	MicroPelletRT <i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Typhimurium</i> WDCM 00121-ATCC® 13311™



© The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalogue marks are trademarks of ATCC. Meccanti S.p.A. is licensed to use these trademarks and sell products derived from ATCC cultures. Look for the ATCC Licensed Derivative Emblem for products derived from ATCC cultures.

# MICROPELLETRTPLUS

code	product name
<b>MPRTP3A017</b>	MicroPelletRTPlus Aspergillus brasiliensis WDCM 00053-ATCC® 16404™, log3
<b>MPRTP3B001</b>	MicroPelletRTPlus Bacillus cereus ATCC® 10876™, log3
<b>MPRTP3B006</b>	MicroPelletRTPlus Geobacillus stearothermophilus ATCC® 7953™, log3
<b>MPRTP3B007</b>	MicroPelletRTPlus Bacillus subsp. Spizizenii WDCM 00003-ATCC® 6633™, log3
<b>MPRTP3B033</b>	MicroPelletRTPlus Bacillus cereus NCIMB 7464, log3
<b>MPRTP3C005</b>	MicroPelletRTPlus Candida albicans WDCM 00054-ATCC® 10231™, log3
<b>MPRTP3C026</b>	MicroPelletRTPlus Clostridium perfringens WDCM 00007-ATCC® 13124™, log3
<b>MPRTP3C029</b>	MicroPelletRTPlus Clostridium sporogenes ATCC® 11437™, log3
<b>MPRTP3C030</b>	MicroPelletRTPlus Clostridium sporogenes WDCM 00008-ATCC® 19404™, log3
<b>MPRTP3E001</b>	MicroPelletRTPlus Enterobacter aerogenes WDCM 00175-ATCC® 13048™, log3
<b>MPRTP3E015</b>	MicroPelletRTPlus Enterococcus faecalis WDCM 00087-ATCC® 29212™, log3
<b>MPRTP3E027</b>	MicroPelletRTPlus Escherichia coli ATCC® 35218™, log3
<b>MPRTP3E032</b>	MicroPelletRTPlus Escherichia coli WDCM 00012,WDCM 00196-ATCC® 8739™, log3
<b>MPRTP3E034</b>	MicroPelletRTPlus Cronobacter muytjensii WDCM 00213- ATCC® 51329™, log3
<b>MPRTP3E041</b>	MicroPelletRTPlus Enterobacter aerogenes NCIMB 10102, log3
<b>MPRTP3E044</b>	MicroPelletRTPlus Escherichia coli ATCC® 51813™, log3
<b>MPRTP3E045</b>	MicroPelletRTPlus Enterococcus faecalis ATCC® 7080™, log3
<b>MPRTP3L013</b>	MicroPelletRTPlus Listeria monocytogenes ATCC® 19115™, log3
<b>MPRTP3L016</b>	MicroPelletRTPlus Lactobacillus fermentum ATCC® 9338™, log3
<b>MPRTP3L021</b>	MicroPelletRTPlus Listeria innocua WDCM 00017-NCTC 11288, log3
<b>MPRTP3L023</b>	MicroPelletRTPlus Listeria monocytogenes WDCM 00021-ATCC® 13932™, log3
<b>MPRTP3M004</b>	MicroPelletRTPlus Kocuria rhizophila ATCC® 9341™, log3
<b>MPRTP3P028</b>	MicroPelletRTPlus Pseudomonas aeruginosa WDCM 00025-ATCC® 27853™, log3
<b>MPRTP3P037</b>	MicroPelletRTPlus Penicillium aurantiogriseum WDCM 00057-ATCC® 16025™, log3
<b>MPRTP3P038</b>	MicroPelletRTPlus Pseudomonas aeruginosa ATCC® 15442™, log3
<b>MPRTP3S018</b>	MicroPelletRTPlus Salmonella enterica subsp. enterica serovar Typhimurium WDCM 00031-ATCC® 14028™, log3
<b>MPRTP3S039</b>	MicroPelletRTPlus Staphylococcus aureus subsp. aureus WDCM 00034-ATCC® 25923™, log3



The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalogue marks are trademarks of ATCC. Mecconti S.A.R.L. Sp. z o. o. is licensed to use these trademarks and sell products derived from ATCC cultures. Look for the ATCC Licensed Derivative Emblem for products derived from ATCC cultures.

# MICROPELLETRTPLUS

code	product name
<b>MPRTP3S047</b>	MicroPelletRTPlus <i>Staphylococcus aureus</i> subsp. <i>aureus</i> WDCM 00032,WDCM 00193-ATCC® 6538™, log3
<b>MPRTP3S048</b>	MicroPelletRTPlus <i>Staphylococcus aureus</i> subsp. <i>aureus</i> WDCM 00033,WDCM 00193,WDCM 00195-ATCC® 6538P™, log3
<b>MPRTP3S051</b>	MicroPelletRTPlus <i>Staphylococcus epidermidis</i> WDCM 00036-ATCC® 12228™, log3
<b>MPRTP3S103</b>	MicroPelletRTPlus <i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Typhimurium</i> WDCM 00121-ATCC® 13311™, log3
<b>MPRTP3S111</b>	MicroPelletRTPlus <i>Staphylococcus aureus</i> subsp. <i>aureus</i> WDCM 00034-NCIMB 12702, log3
<b>MPRTP4B007</b>	MicroPelletRTPlus <i>Bacillus subtilis</i> subsp. <i>Spizizenii</i> WDCM 00003-ATCC® 6633™, log4
<b>MPRTP4E044</b>	MicroPelletRTPlus <i>Escherichia coli</i> ATCC® 51813™, log4
<b>MPRTP4M004</b>	MicroPelletRTPlus <i>Kocuria rhizophila</i> ATCC® 9341™, log4
<b>MPRTP4S018</b>	MicroPelletRTPlus <i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Typhimurium</i> WDCM 00031-ATCC® 14028™, log4
<b>MPRTP4S039</b>	MicroPelletRTPlus <i>Staphylococcus aureus</i> subsp. <i>aureus</i> WDCM 00034-ATCC® 25923™, log4
<b>MPRTP4S051</b>	MicroPelletRTPlus <i>Staphylococcus epidermidis</i> WDCM 00036-ATCC® 12228™, log4
<b>MPRTP4S110</b>	MicroPelletRTPlus <i>Staphylococcus aureus</i> subsp. <i>aureus</i> ATCC® 700699™, log4
<b>MPRTP4S111</b>	MicroPelletRTPlus <i>Staphylococcus aureus</i> subsp. <i>aureus</i> WDCM 00034-NCIMB 12702, log4
<b>MPRTP6A017</b>	MicroPelletRTPlus <i>Aspergillus brasiliensis</i> WDCM 00053-ATCC® 16404™, log6
<b>MPRTP6B007</b>	MicroPelletRTPlus <i>Bacillus subtilis</i> subsp. <i>Spizizenii</i> WDCM 00003-ATCC® 6633™, log6
<b>MPRTP6P028</b>	MicroPelletRTPlus <i>Pseudomonas aeruginosa</i> WDCM 00025-ATCC® 27853™, log6
<b>MPRTP7C005</b>	MicroPelletRTPlus <i>Candida albicans</i> WDCM 00054-ATCC® 10231™, log7
<b>MPRTP7E015</b>	MicroPelletRTPlus <i>Enterococcus faecalis</i> WDCM 00087-ATCC® 29212™, log7
<b>MPRTP7E032</b>	MicroPelletRTPlus <i>Escherichia coli</i> WDCM 00012,WDCM 00196-ATCC® 8739™, log7



\* The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalogue marks are trademarks of ATCC. Meccanti S.p.A. is licensed to use these trademarks and sell products derived from ATCC cultures. Look for the ATCC Licensed Derivative Emblem for products derived from ATCC cultures.

# MICROPELLET RTPLUS

code

product name

**MPRTP7K003** MicroPelletRTPlus Klebsiella pneumoniae subsp. Pneumoniae ATCC® 10031™, log7

**MPRTP7P030** MicroPelletRTPlus Pseudomonas aeruginosa WDCM 00026-ATCC® 9027™, log7

**MPRTP7S047** MicroPelletRTPlus Staphylococcus aureus subsp. aureus WDCM 00032,WDCM 00193-ATCC® 6538™, log7

**MPRTP7S103** MicroPelletRTPlus Salmonella enterica subsp. enterica serovar Typhimurium WDCM 00121-ATCC® 13311™, log7

**MPRTP7S109** MicroPelletRTPlus Salmonella enterica subsp. enterica WDCM 00029-NCTC 6017, log7

**MPRTP8E042** MicroPelletRTPlus Enterobacter cloacae subsp. Cloacae ATCC® 35030™, log8

**MPRTP8S051** MicroPelletRTPlus Staphylococcus epidermidis WDCM 00036-ATCC® 12228™, log8

# MICROPELLET RT10

code

product name

**MPRT10A017** MicroPelletRT10 Aspergillus brasiliensis WDCM 00053-ATCC® 16404™

**MPRT10B007** MicroPelletRT10 Bacillus subtilis subsp. Spizizenii WDCM 00003-ATCC® 6633™

**MPRT10C005** MicroPelletRT10 Candida albicans WDCM 00054-ATCC® 10231™

**MPRT10C026** MicroPelletRT10 Clostridium perfringens WDCM 00007-ATCC® 13124™

**MPRT10C029** MicroPelletRT10 Clostridium sporogenes ATCC® 11437™

**MPRT10C030** MicroPelletRT10 Clostridium sporogenes WDCM 00008-ATCC® 19404™

**MPRT10E032** MicroPelletRT10 Escherichia coli WDCM 00012,WDCM 00196-ATCC® 8739™

**MPRT10P028** MicroPelletRT10 Pseudomonas aeruginosa WDCM 00025-ATCC® 27853™

**MPRT10S018** MicroPelletRT10 Salmonella enterica subsp. enterica serovar Typhimurium WDCM 00031-ATCC® 14028™

**MPRT10S047** MicroPelletRT10 Staphylococcus aureus subsp. aureus WDCM 00032,WDCM 00193-ATCC® 6538™

ATCC Licensed  
Derivative

Zastopnik in distributer: V.I.A. d.o.o. Murska Sobota, Černelavci, Gorička ulica 8, 9000 Murska Sobota,  
Tel.: 02 53 48 406, Fax.: 02 53 48 407, E-mail: [viams@viams.net](mailto:viams@viams.net), [www.viams.net](http://www.viams.net)

**MICROORGANISMS**  
LYOPHILIZED REFERENCE ORGANISMS



MECCONTI S.A.R.L. SP. Z O.O.



**ZASTOPNIK IN DISTRIBUTER**  
**V.I.A. d.o.o. Murska Sobota**  
**Gorička ulica 8, Černelavci**  
**9000 Murska Sobota**  
**Tel: 02 53 48 407**  
**E-mail: info@viams.net**

ZASTOPNIK IN DISTRIBUTER  
V.I.A. d.o.o. Murska Sobota  
Gorička ulica 8, Černelavci  
9000 Murska Sobota  
Tel: 02 53 48 407  
E-mail: [info@viams.net](mailto:info@viams.net)

Zastopnik in distributer V.I.A. d.o.o. Murska Sobota, Černelavci, Gorička ulica 8, 9000 Murska Sobota,  
Tel: 02 53 48 406, Fax: 02 53 48 407, E-mail: [info@viams.net](mailto:info@viams.net), www.viams.net