

# THE BOLT™

## AUTOMATED ANALYZER

### Features & Benefits

- **Modular Design with Optional:**

  - Convection Incubator*

  - Linear Shaker*

  - EIA or EIA+CLIA Reader\**

- Flexible, Open Architecture  
*(Program any EIA or CLIA protocol)*

- Intuitive User Interface

- Space Saving Design  
*(Fits on Standard 60cm lab bench)*

- 96 Sample Capacity

- High Precision Syringes

- Exterior Status Indicator Light

- LIS Connectivity

- Competitive Pricing

- Low Consumable Costs  
*(No disposable tips)*



\*CLIA protocols require EIA+CLIA reader option